

Productivity and robustness

Treats up to 60 samples per hour.
Proven efficiency of blending, including carrot type foods.
2 year manufacturer's guarantee.

Safety

Perfect decontamination: stainless steel inside, access to the whole bowl, dismountable paddles, waste recovery drawer.
Electronic protection of the motor and the blending mechanism.

Quality assurance

Guaranteed quality of blending: optimized microorganism recovery level.

BL 400 Blender

For intensive use without
maintenance




MICROBIOLOGY RANGE

DILUTORS
BLENDERS
MEDIA PREPARATION
MEDIA DISTRIBUTION
BIOCOLLECTORS

Automated Blenders

ACCESSORIES

Bag holder



Sterile irradiated bags

Sterile bags lateral filter

Full page filter bags

Closing clips

Tips Straws 1 ml

Rotary plater



Pulsifier : specially suited to shellfish and meat containing bone fragments

Smasher :
The quietest blender on the market



Dilutor
DI 400

Specifications	BL 400 Blender
Timer	10 sec to 3 mn - continuous
Speed	240 strokes/min
Cleaning	Easy to clean chamber dismountable pads
Waste Drawer	Accessory
Blender Bags	80 to 400 ml
Automatic Start upon closing of the door	Yes
Stainless Steel Interior	Yes
Dimensions WxDxH	260x400x265 mm
Weight	18 kg